

JUN 20 2001

510(k) Summary

| | |
|---------------------|--|
| Introduction | According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence. |
|---------------------|--|

| | |
|--|---|
| 1) Submitter name, address, contact | Roche Diagnostics Corporation 9115 Hague Rd. Indianapolis, IN 46250 (317) 845-2000 Contact Person: Mike Flis Date Prepared: June 4, 2001 |
|--|---|

| | |
|-----------------------|---|
| 2) Device name | Proprietary name: Accu-Chek Active System Classification name: Glucose dehydrogenase, glucose test system (21 C.F.R. § 862.1345)(75LFR) |
|-----------------------|---|

| | |
|----------------------------|---|
| 3) Predicate device | We claim substantial equivalence to the current legally marketed Accu-Chek Simplicity System. |
|----------------------------|---|

| | |
|------------------------------|---|
| 4) Device Description | Instrument Operating Principle -- photometry Reagent Test Principle -- glucose dehydrogenase |
|------------------------------|---|

| | |
|------------------------|--|
| 5) Intended use | The Accu-Chek Active system is designed to quantitatively measure the concentration of glucose in capillary whole blood. The device is indicated for professional use and over-the-counter sale. |
|------------------------|--|

Continued on next page

510(k) Summary, Continued

- 6) **Similarities** The Roche Diagnostics Accu-Chek Active System is substantially equivalent to the current legally marketed Accu-Chek Simplicity System. The proposed modification is relatively modest in scope. The following is a list of some of the claims and features unaffected by the proposed modification.

| Feature/Claim | Detail |
|---|---|
| Intended use | The Accu-Chek Simplicity and Accu-Chek Active systems are designed to quantitatively measure the concentration of glucose in capillary whole blood. The test device is indicated for professional use and over-the-counter sale. This device is not suitable for testing neonate samples. |
| Test principle | Glucose dehydrogenase chemical reaction. The instrument measures the extent of color change (photometric) caused by the presence of glucose in the sample. The amount of color change is related to the glucose concentration in the blood sample. |
| Venous blood samples | Although the method comparison studies provide evidence that the device works acceptably with venous blood samples, the labeling is limited to use with capillary whole blood. |
| Monitor coding procedure | Code chip is provided with each carton of test strips. |
| Test strip storage conditions | Store at room temperature between +36° F(+2° C) and +86° F(+30° C). |
| Test strip operating conditions | Between +5° F(+10° C) and +104° F(+40° C). |
| Quality control testing frequency | Tests should be run with liquid quality control materials whenever a new vial of test strips is opened or an unusual blood test result is obtained. |
| Quality control acceptable range | The mean is strip lot specific and will be determined individually. The range of the controls is within ± 15 mg/dL or $\pm 15\%$ compared to the determined mean. |
| Labeling instructions regarding expected results | The normal fasting adult blood glucose range for a non-diabetic is 70-105 mg/dL. One to two hours after meals, normal blood glucose levels should be less than 140 mg/dL. Doctors will determine the range that is appropriate for the patients. |
| Labeling instructions regarding response to unusual results | Run a quality control test, if the result is outside the acceptable QC recovery range contact Roche Diagnostics's Accu-Chek Customer Care center; if result is within the acceptable range, review proper testing procedure and repeat blood glucose test with a new test strip. |
| Acceptable sample types | Capillary whole blood samples |
| Reportable range | 10-600 mg/dL |
| Hematocrit range | 30 – 55% |
| Warnings and precautions | For <i>in vitro</i> diagnostic use only. |

Continued on next page

510(k) Summary, Continued

**7) Data
demonstrating
substantial
equivalence**

Performance testing on the modified Accu-Chek Active System demonstrated that the device meets the performance requirements for its intended use. A multi-center performance study was conducted to evaluate the accuracy and precision of the modified device. The clinical data demonstrates that the performance of the Accu-Chek Active correlates well with the laboratory plasma glucose reference test method. All predetermined acceptance criteria were satisfied. The data also demonstrates that the Accu-Chek Active is substantially equivalent to the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

JUN 20 2001

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

Mr. Mike Flis
Regulatory Affairs Principal
Roche Diagnostics corporation
9115 Hague Road
Indianapolis, IN 46250

Re: 510(K) Number: K011738
Trade/Device Name: Accu-Chek Active Test System
Regulation Number: 862.1345
Regulatory Class: II
Product Code: NBW, LFR
Dated: June 4, 2001
Received: June 5, 2001

Dear Mr. Flis:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive style with a large, stylized 'S' and 'G'.

Steven I. Gutman, M.D., M.B.A.
Director
Division of Clinical Laboratory Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use Statement

510(k) Number (if known): K011738

Device Name: Accu-Chek Active Test System

Indications for Use:

The Accu-Chek Active system is designed to quantitatively measure the concentration of glucose in capillary whole blood. The device is indicated for professional use and over-the-counter sale.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Fred Loay
(Division Sign-Off)
Division of Clinical Laboratory Devices
510(k) Number K011738

Prescription Use _____
(Per 21 CFR 801.109)

OR

Over-The-Counter Use ☒

(Optional Format 1-2-96)