

MAY 29 2002

K021414

Summary of Safety and Effectiveness  
Liquichek™ Urine Toxicology Control Levels S1 Low Opiate and S2 Low Opiate

1.0 **Submitter**

Bio-Rad Laboratories  
9500 Jeronimo Road,  
Irvine, California 92618-2017  
Telephone: (949) 598-1200  
Fax: (949) 598-1555

**Contact Person**

Ofelia Cachola  
Regulatory Affairs Specialist  
Telephone: (949) 598-1287

**Date of Summary Preparation**

May 1, 2002

2.0 **Device Identification**

Product Trade Name: Liquichek™ Urine Toxicology Control Levels S1 Low Opiate and S2 Low Opiate  
Common Name: Drug Mixture Controls  
Classifications: Class I  
Product Code: DIF  
Regulation Number: 21 CFR 862.3280

3.0 **Device to Which Substantial Equivalence is Claimed**

Liquichek™ Urine Toxicology Control  
Bio-Rad Laboratories  
Irvine, California

Docket Numbers: K992812 and K992809

4.0 **Description of Device**

Liquichek™ Urine Toxicology Control Levels S1 Low Opiate and S2 Low Opiate are prepared from human urine with added constituents of animal origin, drugs, drug metabolites, preservatives, and stabilizers.  
The control is provided in liquid form for convenience.

#### 4. **Statement of Intended Use**

Liquichek™ Urine Toxicology Control is intended for use as quality control urine to monitor the performance of laboratory urine toxicology screening procedures.

#### 5.0 **Comparison of the new device with the Predicate Device**

The new Liquichek™ Urine Toxicology Control claims substantial equivalence to the Liquichek™ Urine Toxicology Control currently in commercial distribution (K992812 and K992809). The new Liquichek™ Urine Toxicology Control has additional claims for Creatinine, pH, and Specific Gravity and the current product does not.

**Table 1.** Similarities and Differences between new and predicate device.

Characteristics	Bio Rad Liquichek™ Urine Toxicology Control (New Device)	Bio Rad Liquichek™ Urine Toxicology Control (Predicate Device)
<b>Similarities</b>		
<b>Intended Use</b>	Liquichek™ Urine Toxicology Control is intended for use as quality control urine to monitor the performance of laboratory urine toxicology screening procedures.	Liquichek™ Urine Toxicology Control is intended for use as a quality control urine to monitor the performance of laboratory urine toxicology screening procedures.
<b>Levels</b>	Same as predicate device.	Level S1 Low Opiate = Drugs added at concentrations 20-25% <b>below</b> immunoassay cutoffs.  Level S2 Low Opiate = Drugs added at concentrations 20-25% <b>above</b> immunoassay cutoffs.
<b>Form</b>	Liquid	Liquid
<b>Matrix</b>	Human urine	Human urine
<b>Storage</b> (Unopened)	2-8° C until expiration date	2-8° C until expiration date
<b>Open Vial Claim</b>	2-8° C for 30 days.	2-8° C for 30 days.
<b>Differences</b>		
<b>Analytes</b>	Same analytes as the predicate device with the additional claims for Creatinine, pH and Specific	D-Amphetamine, Secobarbital, Nordiazepam, 11-Nor-Δ-9-THC-9-COOH), Benzoyllecgonine,

	Gravity.	Ethanol, Lysergic Acid Diethylamide (LSD), Methadone, Methaqualone, Morphine (Free), Phencyclidine, Propoxyphene, Nortriptyline.
--	----------	--

## 2.0 STATEMENT OF SUPPORTING DATA

Stability studies have been performed to determine the open vial stability and shelf life for the Liquichek™ Urine Toxicology Control. Product claims are as follows:

2.1 Open vial: Once the control is opened, all analytes will be stable for 30 days when stored tightly capped at 2-8°C.

2.2 Shelf Life: 24 months when stored at 2-8°C

Real time studies will be ongoing to support the shelf life of this product.

All supporting data is retained on file at Bio-Rad Laboratories.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

**MAY 29 2002**

Ms. Elizabeth Platt  
Regulatory Affairs Manager  
Bio-Rad Laboratories, QSD  
9500 Jeronimo Road  
Irvine, California 92618-2017

Re: k021414  
Trade/Device Name: Liquichek™ Urine Toxicology Control Levels S1 Low Opiate and  
S2 Low Opiate  
Regulation Number: 21 CFR § 862.3280  
Regulation Name: Clinical Toxicology Control Material  
Regulatory Class: I  
Product Code: DIF  
Dated: May 1, 2002  
Received: May 3, 2002

Dear Ms. Platt:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive, slightly slanted style.

Steven I. Gutman, M.D., M.B.A.  
Director  
Division of Clinical  
Laboratory Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health


Enclosure

510 (k) Number (if known): K021414

Device Name: **Liquichek™ Urine Toxicology Control Levels S1  
Low Opiate and S2 Low Opiate**

Indications for Use:

**A quality control urine to monitor the performance of laboratory urine  
toxicology screening procedures.**

  
(Division Sign-Off)  
Division of Clinical Laboratory Devices  
510(k) Number K021414

(PLEASE DO NOT WRITE BELOW THE LINE-CONTINUE ON ANOTHER PAGE IF  
NEEDED)

---

*Concurrence of CDRH, Office of Device Evaluation (ODE)*

Prescription use X or Over-the Counter  
use \_\_\_\_\_