

K011621

**AUG 1 0 2001**

**510(k) Summary**

**Substantial Equivalence & Safety & Effectiveness Information**

[May 18, 2001]

- I. Proprietary Name:** Koala Clamp & Cutter
- II. Common Name:** Umbilical Cord Clamp & Cutter Device
- III. Classification Name:** Umbilical Clamp (HFW), Umbilical Scissors (HDJ)
- IV. Substantial Equivalence:** The Koala Clamp & Cutter is substantially equivalent to the ClampCut by PriceInvena APS, a legally marketed device in the US. [510(k) No. 982464]
- V. Statement of Intended Use:** The indications for use of the Koala Clamp & are to simultaneously clamp and cut the umbilical cord.

KO 11621

**VI. Comparison of Present Device & Predicate Device:**

<b>Comparison Chart of Koala Clamp &amp; Cutter and Predicate Device</b>		
	<b><i>Koala Clamp</i></b>	<b><i>ClampCut (SCC-23)</i></b>
<b><i>Indications for Use</i></b>	The indications for use are to simultaneously clamp and cut the umbilical cord.	The indications for use are to simultaneously clamp and cut the umbilical cord. The device reduces the risk to healthcare practitioners of unnecessary exposure to infection by bloodborne diseases.
<b><i>Target Population</i></b>	All births	All births
<b><i>Materials: Plastic</i></b>	Lexan 144R (alternate: Lexan 124R)	Polyamid 6.6.
<b><i>Biocompatible</i></b>	Yes	Yes
<b><i>Materials: Blade</i></b>	Stainless Steel Surgical Blade	Stainless Steel Knife
<b><i>Shield (for Blade)</i></b>	Yes	No
<b><i>Dimensions: Open Unit</i></b>	74mm x 50.6mm x 32.9mm	120mm x 70mm x 20mm
<b><i>Dimensions: Closed Unit</i></b>	46.8mm x 50.6mm x 32.9mm	85mm x 60mm x 20mm
<b><i>Dimensions: Infant Clamp</i></b>	43.5mm x 23.9mm x 6.5mm	50mm x 9mm x 7mm
<b><i>Weight of Total Unit</i></b>	1.316 oz	0.811 oz (23 grammes)
<b><i>Weight of Cutter Unit alone</i></b>	0.53 oz	0.713 oz (20.2 grammes)
<b><i>Weight of Infant Clamp alone</i></b>	0.086 oz	0.099 oz (2.8 grammes)
<b><i>Sterility</i></b>	Gamma-Sterilized (alternate: Ethylene Oxide Gas)	Ethylene Oxide Gas
<b><i>Anatomical Sites</i></b>	Umbilical Cord	Umbilical Cord



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

AUG 1 0 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Ronald B. Hicks  
Vice President and COO  
Maternus, Inc.  
P.O. Box 782089  
SAN ANTONIO TX 78278-2089

Re: K011621  
Koala Clamp & Cutter, Umbilical Cord Clamp and Cutter  
Dated: May 22, 2001  
Received: May 25, 2001  
Regulatory Class: II  
21 CFR §884.4530/Procode: 85 HFW

Dear Mr. Hicks:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4639. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure (s)

## Statement of Indications for Use


510(k) Number: K011621


Device Name: Koala Clamp & Cutter

Indications for Use: The indications for use of the Koala Clamp & Cutter are to simultaneously clamp and cut the umbilical cord.

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use   
(Per 21 CFR 801.109)

  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number K011621