


K011749

AUG 22 2001

Siam Sempermed Corp., Ltd. 
110 Moo 8 Kanjanavanit Rd., Hat Yai, Songkhla, Thailand 90230
Tel: 66 074 291 648 to 9 Fax: 66 074 291 650

Page 1 (2)

510 (k) SUMMARY
21062001

1.0 APPLICANT:

Dr. POONSUK CHERDKIATGUMCHAI
SIAM SEMPERMED CORPORATION., Ltd.
110 MOO 8 KANJANAVANIT ROAD
PATHONG HATYAI SONGKHLA
THAILAND 90230
TEL: 66 074 291 648 OR 291 649
FAX: 66 074 291 650

2.0 CONTACT PERSON

Dr. POONSUK CHERDKIATGUMCHAI
SIAM SEMPERMED CORPORATION., Ltd.
110 MOO 8 KANJANAVANIT ROAD
PATHONG HATYAI SONGKHLA
THAILAND 90230
TEL: 66 074 291 648 OR 291 649
FAX: 66 074 291 650

MR WILLIAM HARRIS
SEMPERMED USA Inc.
30798 US Hwy. 19 N
Palm Harbor,
USA FL 34684
TEL: 727 787 7250
FAX: 727 787 7558

3.0 Device Class: I

Product code: 80LYY

4.0 Specification: Latex patient examination glove, Powder Glove (Single side polymer coated) -Class I 80LYY
meets all of the requirements of ASTM standard D3578-00

5.0 Device Description: Latex Patient Examination glove, Powder Glove (Single side polymer coated), non sterile
200 micrograms or less of total water extractable protein per gram


6.0 Intended use: A patient examination glove is a disposable device intended for medical purposes that is worn
on the examiners hand or finger to prevent contamination between patient and examiner.

7.0 Outer Surface : Free from talc (Magnesium silicate)

8.0 Primary Dermal Irritation in Rabbits Guinea Pig Sensitization (Buehler) : Consumer Product Testing Co.
Experiment reference number : T95-0189-1

Conclusion : According to Federal Hazardous Substances Act Regulation, (16 CFR 1500.41), and under the
conditions of this test, This test article is not a primary dermal irritant
: This test article is not a sensitizer in guinea-pigs, under condition of this test.

K011749

Siam Sempermed 
Corp., Ltd.

Page 2 (2)

510 (k) SUMMARY
21062001

9.0 QUALITY CHARACTERISTICS

Dimensions	Meet ASTM D 3578-00
Physical Properties	Meet ASTM D 3578-00
Protein content	Recommended 200 ug/dm ² in ASTM D 3578-00
Freedom from pinholes	Meet ASTM D 3578-00 Meet ASTM D 5151

10. Conclusion: Siam Sempermed Latex Patient Examination Glove ,Powder Glove , 200 micrograms or less of total water extractable protein per gram
meet the ASTM standard or equivalent standard
meet pinhole FDA requirements
meet labeling claims (see 5.0 and 6.0 above)

P. Cherdkiatgumchai

Dr. POONSUK CHERDKIATGUMCHAI
Chief Quality Officer
21062001



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

AUG 22 2001

Ms. Katie Levinson
Product Manager
Sempermed USA, Incorporated
30798 US Highway 19 North
Palm Harbor, Florida 34684

Re: K011749
Trade/Device Name: Latex Powdered Patient Examination
Glove, 200 Micrograms or Less
Regulation Number: 880.6250
Regulatory Class: I
Product Code: LYY
Dated: May 24, 2001
Received: May 31, 2001

Dear Ms. Levinson:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements

concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, appearing to read "Timothy A. Ulatowski" with a stylized flourish at the end.

Timothy A. Ulatowski
Director
Division of Dental, Infection Control
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

INDICATION FOR USE

Applicant: Siam Sempermed Corp. Ltd.

510(k) Number (if known): K011749


Device Name: Latex Powdered Examination Glove with a Protein Content Labeling
Claim of 200 micrograms or less of water extractable protein per glove —

Indications For Use:

A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hands or finger to prevent contamination between patient and examiner.
(21CFR 880.6250)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)

Division of Dental, Infection Control,
and General Hospital Devices

510(k) Number K011749