

NOV 05 2001

F. 510(k) Summary

Applicant: Stick Tech Ltd, PO Box 114, 20521 Turku, Finland
Tel: +358-2-4808 2500
Fax: +358-2-241 0032

Contact Person: Ilkka Kangasniemi, Ph.D.

U.S. Agent to respond to
FDA requests: William M. Troetel, Ph.D.
80 Parkway West
Mount Vernon, NY 10552
Tel: (914) 664-1640
Fax: (914) 667-3993
wtroetel@westnet.com

Date Prepared: June 4th, 2001

Device Trade Name: everStickNet™

Device Common Name : Glass fiber reinforcement material

Device Classification Name: Denture relining, repairing, or rebasing resin
(21 CFR §872.3760)

Description of Device:

everStickNet™ is a semi-manufactured product made of glass fibers and polymer/resin matrix for reinforcing dental acrylic polymers. everStickNet™ is made of a thin fiberglass fabric, which increase the strength and stiffness of the final product in all directions.

Intended Use:

As reinforcement in manufacturing and/or repairing full or partial dentures as well as overdentures and orthodontic appliances.

As reinforcement for temporary and/or permanent plastic/composite partial and full crowns and bridges.

As reinforcement for customized splints used to immobilize teeth which may be required post-trauma, post-operative, or for orthodontic therapy.

everStickNet™ is substantially equivalent to fibreNet™, approved under 510(k) number K003334 dated January 8, 2001.

Testing which has been performed on everStickNet™ indicates that the device has the same intended use but somewhat different technological characteristics.

everStickNet™ is a polymer impregnated bi-directionally continuously woven glass fiber whereas fibreNet™ is a polymer pre-impregnated bi-directionally woven glass fiber. The different technological characteristics of everStickNet™ does not raise new questions of safety and effectiveness and demonstrates that the device is as safe and effective as the predicate device.

Test results indicate that there are no hazards presented with the use of everStickNet™ as compared with the predicate device.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Stick Tech Limited
C/O Dr. William Troetel
Managing Partner
William M. Troetel, LLC
80 Parkway West
Mount Vernon, New York 10552

NOV 05 2001

Re: K011799

Trade/Device Name: Eversticknet
Regulation Number: 872.3760
Regulation Name: EBI Resin, Denture, Relining, Repairing, Rebasing
Regulatory Class: II
Product Code: EBI
Dated: October 29, 2001
Received: October 30, 2001

Dear Dr. Troetel:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements


of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4613. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Timothy A. Ulatowski
Director
Division of Dental, Infection Control
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

K011799

B-1

NOV 05 2001


B. Indications for Use

510(k) Number (if known): K011799

Device Name: everStickNet™

Indications for Use:

- As reinforcement in manufacturing and/or repairing full or partial dentures as well as overdentures and orthodontic appliances.
- As reinforcement for temporary and/or permanent plastic/composite partial and full crowns and bridges.
- As reinforcement for customized splints used to immobilize teeth which may be required for post-trauma, post-operative, or orthodontic therapy.



(Division Sign-Off)

Division of Dental, Infection Control,
and General Hospital Devices

510(k) Number K011799

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801.109)

OR

Over-the-Counter Use

(Optional Format 1-2-96)