

AUG 22 2001

Summary of Safety and Effectiveness

K011856

Encore Orthopedics[®], Inc.
9800 Metric Blvd
Austin, TX 78758
512-832-9500

Trade Name: Radiographic Markers

Classification Name: Marker, radiographic, implantable

Description: The Radiographic Markers are spherical in shape and available in Ø0.8-1.0 mm. These markers are manufactured from tantalum conforming to the composition called out in ASTM F560-98.

Intended Use: The Radiographic Markers are intended for implantation in the bone to aid the surgeon in radiographic assessment of total hip or knee component movement.

Comparable Features to Predicate Device(s): Features comparable to predicate devices include same materials, design and indications as the Biomet Tantalum Beads (K010348).



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

AUG 22 2001

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Ms. Joanne Droege
Regulatory/Quality Assurance Engineer
Encore Orthopedics, Inc.
9800 Metric Boulevard
Austin, Texas 78758

Re: K011856
Trade/Device Name: Radiographic Markers
Regulation Number: 878.4300
Regulatory Class: II
Product Code: NEU, FZP
Dated: June 11, 2001
Received: June 13, 2001

Dear Ms. Droege:


We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, appearing to read "for Mark N. Melanson", is written over the typed name of the signatory.

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative
and Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

K011856

510(k) Number (if known): _____

Device Name: Radiographic Markers

Indications For Use:

Radiographic Marker
Indications For Use

The Radiographic Markers are intended for implantation in the bone to aid the surgeon in radiographic assessment of total hip or knee component movement. The Radiographic Markers are manufactured from inert, biocompatible materials.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(per 21 CFR 801.109)

OR Over-The-Counter Use _____

(Optional Format 1-2-96)_____

for Mark A. Miller
(Division Sign-Off)

Division of General, Restorative
and Neurological Devices

510(k) Number K011856

RA
II

SK23