

JUL - 9 2001

K011864  
p1 of 3

## SUMMARY OF SAFETY AND EFFECTIVENESS

### COMPANY AND CONTACT PERSON

Medtronic Perfusion Systems  
7611 Northland Drive  
Minneapolis, MN 55428  
Tel: (763) 391-9136  
FAX: (763) 391-9603

Stephani K. Ayala, Regulatory Affairs Specialist, Regulatory Affairs

### DEVICE NAME

Trillium™ MYOthem™ XP Cardioplegia Delivery System

### NAME OF PREDICATED OR LEGALLY MARKETING DEVICE

MYOthem™ XP Cardioplegia Delivery System (K971105)  
AFFINITY® Hollow Fiber Oxygenator with Trillium™ Biopassive Surface (K973760)

### DESCRIPTION OF DEVICE

The Trillium™ MYOthem™ XP Cardioplegia Delivery System is designed to mix arterial blood from oxygenators with asanguineous cardioplegia solution. The blood/cardioplegia solution is then cooled/warmed and delivered to the patient.

### STATEMENT OF INTENDED USE

The Trillium™ MYOthem™ XP Cardioplegia Delivery System is a device intended for the mixing, warming/cooling and delivery of oxygenated blood/cardioplegia solution in a predetermined ratio. Blood/cardioplegia solution is delivered to the patient through the cardioplegia delivery system and appropriate cannula by the operation of a single occlusive roller pump.

### STATEMENT OF INTENDED USE OF PREDICATED/MARKETING DEVICE

The MYOthem™ XP Cardioplegia Delivery System is a device intended for the mixing, warming/cooling and delivery of oxygenated blood/cardioplegia solution in a predetermined ratio. Blood/cardioplegia solution is delivered to the patient through the cardioplegia delivery system and appropriate cannula by the operation of a single occlusive roller pump.

### STATEMENT OF TECHNOLOGICAL CHARACTERISTICS COMPARISON

Information regarding technological characteristics comparison is provided in the following section, "Determination of Substantial Equivalence".

**DETERMINATION OF SUBSTANTIAL EQUIVALENCE**

This "SPECIAL 510(k)" is being submitted for a modification to the MYOthem™ XP Cardioplegia System. The modification to the current MYOthem™ XP Cardioplegia System is to coat the blood contact surfaces with Trillium™.

The Trillium™ MYOthem™ XP Cardioplegia Delivery System is being compared to the following Marketed Devices:

- MYOthem™ XP Cardioplegia System (K971105)
- AFFINITY® Hollow Fiber Oxygenator with Trillium™ Biopassive Surface (K973760)

The Trillium™ MYOthem™ XP Cardioplegia Delivery System has the same indications statement and intended uses as the:

- MYOthem™ XP Cardioplegia System (K971105)

The Trillium™ MYOthem™ XP Cardioplegia Delivery System has "no new technological characteristics (e.g., materials and manufacturing processes)" from the MYOthem™ Cardioplegia System. The technological characteristic is solely the coating material of the blood pathway:

- Trillium™

The technological characteristic of the Trillium™ Biopassive Surface is common to other hollow fiber oxygenators currently in commercial distribution as follows:

- AFFINITY® Hollow Fiber Oxygenator with Trillium™ Biopassive Surface (K973760)

This technological characteristic "could affect the safety and effectiveness of the device". However, these "new technological characteristics do not raise new types of safety or effectiveness questions". In addition, "there are acceptable scientific methods which exist for assessing effects of these new technological characteristics".

"Performance data to assess the effects of these new technological characteristics" has been performed. These "performance data demonstrate" that the Trillium™ MYOthem™ XP Cardioplegia Delivery System is substantially equivalent to other marketed cardioplegia delivery systems.

The biocompatibility and *in vitro* bench testing demonstrated that when compared to the predicate devices, the Trillium™ MYOthem™ XP Cardioplegia Delivery System does not significantly affect safety and effectiveness and are substantially equivalent to other commercially distributed cardioplegia delivery system. The *in vitro* bench testing included analysis of:

#### Coating Characteristics

- Coating
- Leaching

#### Physical Characteristics

- Pressure integrity
- Priming volume

#### Performance Characteristics

- Ease of prime
- Heat exchanger performance
- Pressure drop
- Blood trauma



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

JUL - 9 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Medtronic Perfusion Systems  
c/o Ms Marie Holm  
Associate Product Regulations Manager  
7611 Northland Drive N  
Minneapolis, MN 55428.1088

Re: K011864

Trade Name: Trillium™ MYOotherm™ XP Cardioplegia Delivery System (Models: XP41T  
and XP41BT)

Regulation Number: 870.4400

Regulatory Class: II (Two)

Product Code: DTN

Dated: July 5, 2001

Received: July 5, 2001

Dear Ms. Holm:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might


Page 2 - Ms Marie Holm

have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4648. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, appearing to read 'James E. Dillard III', written in a cursive style.

James E. Dillard III  
Director  
Division of Cardiovascular and  
Respiratory Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

**Indications for Use**

**K011864**

**Device Name: Trillium™ MYOtherm™ XP Cardioplegia Delivery System**

**Indications for Use:**

The Trillium™ MYOtherm™ XP Cardioplegia Delivery System is a device intended for the mixing, warming/cooling and delivery of oxygenated blood/cardioplegia solution in a predetermined ratio. Blood/cardioplegia solution is delivered to the patient through the cardioplegia delivery system and appropriate cannula by the operation of a single occlusive roller pump.


---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒  
Per 21 CFR 801.109

OR

Over-The-Counter-Use

  
Division of Cardiovascular & Respiratory Devices  
510(k) Number K011864