

AUG 29 2001

K011 984

I. 510(k) Summary of Safety and Effectiveness

510(k) Summary Of Safety and Effectiveness

I. General Information

This Summary of Safety and Effectiveness information is being submitted in accordance with the requirements of the SMDA of 1990 and 21 § 807.92

Establishment:

- Address: Becton Dickinson VACUTAINER Systems
1 Becton Drive
Franklin Lakes, NJ 07417-1885
- Registration Number: 2243072
- Contact Person: M. Wendy Bosshardt
Regulatory Affairs Specialist
Telephone no.: 201-847-6280
Fax No. 201-847-4858
- Date of Summary: June 25, 2001

Device

- Trade Name: BD VACUTAINER™ Push Button Blood Collection Set
- Classification Name: Tubes, Vials, Systems, Serum Separators, Blood Collection
- Classification: Class II
- Performance Standards: None Established under 514 of the Food, Drug and Cosmetic Act

II. Safety and Effectiveness Information Supporting Substantial Equivalence

Substantial Equivalence Declaration:

The term "Substantial Equivalence" as used in this 510(k) Premarket Notification is limited to the definition of Substantial Equivalence found in the Federal Food, Drug, and Cosmetic Act, as amended and as applied under 21 CFR § 807, Subpart E, under which a device can be marketed without pre-market approval or reclassification. A determination of substantial equivalency under this notification is not intended to have any bearing whatsoever on the resolution of patent infringement suits or any other patent matters. No statements related to, or in support of, substantial equivalence herein shall be construed as an admission against interest under the US Patent Laws or their application by the courts.

- **Device Description**

The BD VACUTAINER™ Push Button Blood Collection Set is for venous blood collection. The wing set contains a needle that will retract into the body of the device when a button is depressed, helping to prevent accidental needle sticks. The retraction of the needle occurs when the user depresses the button.

- **Intended Use**

The BD VACUTAINER™ Push Button Blood Collection Set is a sterile, multiple-sample, single-use winged blood collection set intended for venipuncture to obtain blood specimens from patients. The recommended use of the device is to activate the needle prior to removal from the venipuncture site. The retraction of the intravenous (IV) end of the needle aids in the prevention of needlestick injury.

- **Synopsis of Performance Study Results**

Mechanical and Simulated Use testing was performed to confirm the robustness of the design and reliability of the safety feature function. Clinical use testing was performed to evaluate the function of the safety feature in a blood drawing environment.

III. Predicate Device Summary Table

- Substantial Equivalence

Based on comparison of the device features, materials, intended use and performance, the BD VACUTAINER™ Push Button Blood Collection Set be shown to be substantially equivalent to the commercially available predicate devices indicated in the table below. The predicate devices, K number, and clearance date are also identified in the table below.

Manufacturer	Predicate Device	K-Number	Clearance Date
Becton Dickinson	BD VACUTAINER™	K921636	April 29, 1992
	Brand Safety-Lok™	K931367	May 18, 1993
	Blood Collection Set	K965202	February 10, 1997
		K980414	March 3, 1998

M. Wendy Bosshardt
M. Wendy Bosshardt
Regulatory Affairs Specialist
Becton Dickinson VACUTAINER Systems
Becton Dickinson and Company

August 17, 2001
Date



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

AUG 29 2001

Ms. M. Wendy Bosshardt
Regulatory Affairs Specialist
Becton Dickinson & Company
1 Becton Drive
Franklin Lakes, New Jersey 07417-1880

Re: K011984
Trade/Device Name: BD Vacutainer™ Push Button
Blood Collection Set
Regulation Number: 880.5570 and 862.1675
Regulatory Class: II
Product Code: FMI and JKA
Dated: June 25, 2001
Received: June 26, 2001

Dear Ms. Bosshardt:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements

concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4618. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski
Director
Division of Dental, Infection Control
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

B. INDICATIONS FOR USE

510(k) NUMBER (IF KNOWN): K 011984

DEVICE NAME: BD VACUTAINER™ PUSH BUTTON BLOOD COLLECTION SET

INDICATIONS FOR USE:

The BD VACUTAINER™ Push Button Blood Collection Set is a sterile, multiple-sample, single-use winged blood collection set intended for venipuncture to obtain blood specimens from patients. The recommended use of the device is to activate the needle prior to removal from the venipuncture site. The retraction of the intravenous (IV) end of the needle aids in the prevention of needlestick injury.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)


CONCURRENCE OF CDRH, OFFICE OF DEVICE EVALUATION (ODE)

PRESCRIPTION USE _____

OR OVER-THE-COUNTER USE _____

(PER 21 CFR § 801.109)

(OPTIONAL FORMAT 1-2-96)


(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices
510(k) Number K 011984