

JUL 26 2001

K012095

Implex Corp.

LPS - Continuum<sup>®</sup> TM Hybrid Tibia  
Special 510(k) Premarket Notification

## **510(k) SUMMARY OF SAFETY AND EFFECTIVENESS**

### **The LPS - Continuum TM Hybrid Tibia**

**Submitter Name:** Implex Corp.

**Submitter Address:** 80 Commerce Drive  
Allendale, New Jersey 07401-1600

**Contact Person:** Robert Poggie

**Phone Number:** (201) 818-1800

**Fax Number:** (201) 995-9722

**Date Prepared:** June 29, 2001

**Device Trade Name:** The LPS – Continuum TM Hybrid Tibia

**Device Common Name:** Articular Surface and Tibial Components

**Classification Number and Name:** 21 CFR § 888.3560

---

**Substantial Equivalence:** The term “substantial equivalence” as used in this 510(k) notification is limited to the definition of substantial equivalence found in the Federal Food, Drug and Cosmetic Act, as amended and as applied under 21 CFR 807, Subpart E under which a device can be marketed without premarket approval or reclassification. A determination of substantial equivalency under this notification is not intended to have any bearing whatsoever on the resolution of patent infringement suits or any other patent matters. No statements related to, or in support of substantial equivalence herein shall be construed as an admission against interest under the US Patent Laws or their application by the courts.

**Device Description:** The LPS – Continuum TM Hybrid Tibia is manufactured from Trabecular Metal (Hedrocel Porous Tantalum) with direct compression molded ultra-high molecular weight polyethylene (UHMWPE).

These tibial components are intended for use with Zimmer NexGen LPS Femoral components.

### ***510(k) Summary (Continued)***

**Indications for Use:** The LPS - Continuum<sup>®</sup> TM Hybrid Tibia is intended for use where severe degeneration, trauma, or other pathology of the knee joint indicates cemented total knee arthroplasty. In addition, the LPS - Continuum<sup>®</sup> TM Hybrid Tibia is indicated for use in the presence of knee instability caused by a compromised or non-functional posterior cruciate ligament.

**Conclusion:** The LPS - Continuum TM Hybrid Tibia is substantially equivalent to the identified predicate devices.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

**JUL 26 2001**

Dr. Robert A. Poggie  
Director of Applied Research  
Implex Corp.  
80 Commerce Drive.  
Allendale, New Jersey 07401

Re: K012095  
Trade Name: The LPS-Continuum™ Hybrid Tibia  
Regulatory Class: II  
Regulation Number: 888.3560  
Product Code: JWH  
Dated: July 3, 2001  
Received: July 5, 2001

Dear Dr. Poggie:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general control provisions of the Act. The general control provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, appearing to read 'C. Witten', is written over the typed name.

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative  
and Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

510(k) Number (if  
known):K012095

Device Name:

The LPS- Continuum<sup>®</sup> TM Hybrid Tibia

Indications For Use:

The LPS - Continuum<sup>®</sup> TM Hybrid Tibia is intended for use where severe degeneration, trauma, or other pathology of the knee joint indicates cemented total knee arthroplasty. In addition, the LPS - Continuum<sup>®</sup> TM Hybrid Tibia is indicated for use in the presence of knee instability caused by a compromised or non-functional posterior cruciate ligament.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF  
NEEDED)

Concurrence of CDRH; Office of Device Evaluation (ODE)

Prescription  
Use ✓  
(Per 21 CFR 801.109)

Robert J. Lee MD for CDRH  
(Division Sign-Off)  
Division of General, Restorative  
and Neurological Devices

Counter Use  
510(k) Number K012095

(Optional Format 1-2-96)