

OCT 11 2001

510(k) Summary of Safety and Effectiveness

K012315

Trade Name	LeVeen™ Electrode
Common/Classification Name	Electrosurgical Electrode
Classification	Class II
Summary of Substantial Equivalence	The modified LeVeen™ Electrode is substantially equivalent to the previously cleared LeVeen™ Electrode
Device Description	The device consists of a preshaped, multi-armed electrode array which is contained within a delivery cannula. The array is attached to a handle mechanism that deploys the array into the targeted tissue. In addition, each device is supplied with an insulated introducer sheath/stylet assembly to allow the physician to locate the target tissue prior to placing and deploying the electrode. The device is connected to a RadioTherapeutics RF generator so that RF energy passes from the array to a patient ground pad and heats the tissue surrounding the array.
Intended Use	The LeVeen™ Electrode is intended to be used in conjunction with a RadioTherapeutics Corporation radiofrequency (RF) generator for the thermal coagulation necrosis of soft tissue, including partial or complete ablation of nonresectable liver lesions.

July 20, 2001

RadioTherapeutics Corporation
1308 Borregas Avenue, Sunnyvale, CA 94089
(408) 745-3200 • Fax (408) 745-9848



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

OCT 11 2001

Ms. Kirsten Valley
Vice President, Regulatory
RadioTherapeutics Incorporated
1308 Borregas Avenue
Sunnyvale, California 94089

Re: K012315

Trade/Device Name: LeVeen™ Electrode
Regulation Number: 878.4400
Regulation Name: Electrosurgical cutting and coagulation device and accessories
Regulatory Class: II
Product Code: GEI, JOS
Dated: July 20, 2001
Received: July 23, 2001

Dear Ms. Valley :

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): ~~Unknown~~

K 012315

Device Name: The LeVeen™ Electrode

Indications for Use: The LeVeen™ Electrode is intended to be used in conjunction with a RadioTherapeutics Corporation radiofrequency (RF) generator for the thermal coagulation necrosis of soft tissue, including partial or complete ablation of nonresectable liver lesions.

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)
Division of General, Restorative
and Neurological Devices

510(k) Number K012315

Prescription Use ☒

OR

Over-The-Counter Use ☐

(Optional Format 1-2-96)