

JAN 22 2002

K012480

## 510(k) Summary of Safety & Effectiveness

---

<b>Submitter</b>	Vanguard Medical Concepts, Inc. 5307 Great Oak Drive Lakeland, FL 33815
------------------	---

---

<b>Contact</b>	Mr. Mike Sammon, Ph.D. Director, Research and Development (863) 683-8680, extension 228 (801) 327-3339 (facsimile) mikes@safe-reuse.com
----------------	---

---

<b>Date</b>	July 27, 2001
-------------	---------------

---

<b>Device</b>	<ul style="list-style-type: none"><li>• Trade Names: Vanguard Reprocessed Inflation Devices<ul style="list-style-type: none"><li>⇒ ACS Indeflator™ Inflation Devices</li><li>⇒ Medtronic Everest™ 20/30 Inflation Devices</li><li>⇒ Merit Medical Systems Monarch™/ IntelliSystem®25 Inflation Devices</li><li>⇒ SCIMED® Encore® 26 Inflation Devices</li></ul></li><li>• Common Name: Inflation device or syringe</li><li>• Classification: Class II – Syringe, Balloon Inflation</li><li>• Product Code MAV</li></ul>
---------------	---

---

<b>Predicate Devices</b>	Respective ACS, Medtronic, Merit Medical, and SCIMED® legally marketed inflation devices under various 510(k) premarket notifications.
--------------------------	--

---

<b>Indications for Use</b>	The inflation device is designed to inflate and deflate an angioplasty balloon or other interventional device and to monitor pressure within the balloon.
----------------------------	---

---

*Continued on next page*

## 510(k) Summary of Safety & Effectiveness, Continued

---

<b>Device Description</b>	<p>An inflation device is designed to generate and monitor pressure and is used in angioplasty or other interventional procedures. It is a sterile hand-held mechanical device comprised of a 20 cc syringe, pressure gauge, plunger, connector tube and a male luer fitting for connection to a catheter. The device may be used in conjunction with other accessories such as a hemostatic valve or stopcock, guide wire introducer or torque device.</p> <p>Vanguard receives previously used inflation devices from healthcare facilities; cleans, inspects, tests, repackages and sterilizes the devices; and returns them to the healthcare facility.</p>
<b>Technological Characteristics</b>	<p>The Vanguard reprocessed inflation devices are essentially identical to the currently marketed OEM devices. No changes are made to the currently marketed device's specifications and they possess the same technological characteristics. Biocompatibility and performance/functional testing demonstrate that the devices are equivalent and continue to be safe and effective for their intended use.</p>
<b>Test Data</b>	<p>Sterilization and packaging validations, and functional/performance, and biocompatibility testing demonstrates that the reprocessed devices perform as intended and are safe and effective.</p>
<b>Conclusion</b>	<p>Based on the information provided herein and the 510(k) "Substantial Equivalence" Decision Making Process Chart, we conclude that the Vanguard reprocessed inflation devices are substantially equivalent to the predicate devices under the Federal Food, Drug and Cosmetic Act.</p>

---



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JAN 22 2002

Mike Sammon, Ph.D.  
Director, Research and Development  
Vanguard Medical Concepts, Inc.  
5307 Great Oak Drive  
Lakeland, FL 33815

Re: K012480  
Vanguard Reprocessed Inflation Devices  
Regulation Number: 870.1650  
Regulation Name: Angiographic Injector and Syringe  
Regulatory Class: Class II  
Product Code: 74 MAV  
Dated: November 25, 2001  
Received: November 27, 2001

Dear Dr. Sammon:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with


Page 2 - Mike Sammon, Ph.D.

all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4586. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

  
Bram Zuckerman, M.D.  
Acting Director  
Division of Cardiovascular and  
Respiratory Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number: K012480

Device Name: Vanguard Reprocessed Inflation Devices

### Indications for Use:

The inflation device is designed to inflate and deflate an angioplasty balloon or other interventional device and to monitor pressure within the balloon.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED.)


-----  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use ☐

(Optional Format 1-2-96)

  
Division of Cardiovascular & Respiratory Devices  
510(k) Number K012480