

NOV 1 9 2001

510(k) SUMMARY

K012806

Submitter: Cynsoure, Inc.
10 Elizabeth Drive
Chelmsford, MA 01824

Contact: George Cho
Senior Vice President of Medical Technology

Date Summary Prepared: August 17, 2001

Device Trade Name: PhotoGenica VL

Common Name: Medical Laser System

Classification Name: Instrument, surgical, powered, laser
79-GEX
21 CFR 878.48

Equivalent Device: NLite Laser by ICN Pharmaceuticals, Inc.

Device Description: The PhotoGenica VL is a pulse-dye laser, having the organic dye as the lasing medium. It is a pulsed laser with a wavelength of 580 to 590nm.

Laser activation is both by finger switch and footswitch. Overall weight of the laser is 285lbs, and the size is 44"x19"x24" (HxWxD).

Electrical requirement is 110 VAC or 220 VAC, 20A, 50-60 Hz, single phase.

Intended Use: The PhotoGenica VL is indicated for treatment of periocular wrinkles and dermatological lesions.

Comparison: The PhotoGenica VL laser has an equivalent indication for uses, the same principle of operation, the same wavelength and pulse energy range as the predicate devices.

Nonclinical Performance Data: none

Clinical Performance Data: none

Conclusion: The PhotoGenica VL Laser is another safe and effective device for dermatologic applications.

Additional Information: none



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. George Cho
Senior Vice President
of Medical Technology
Cynosure, Inc.
10 Elizabeth Drive
Chelmsford, Massachusetts 01824

NOV 19 2001

Re: K012806
Trade/Device Name: PhotoGenica VL Dermatology Laser
Regulation Number: 878.4810
Regulation Name: Laser surgical instrument for use in general and
plastic surgery and in dermatology.
Regulatory Class: II
Product Code: GEX
Dated: August 17, 2001
Received: August 21, 2001

Dear Mr. Cho:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050. .

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

501(k) Number (if known): K 012806

NOV 19 2001

Device Name: Cynosure PhotoGenica VL

Indications For Use:

The Cynosure PhotoGenica VL laser is indicated for use in Dermatological and Plastic Surgery applications and the treatment of periocular wrinkles.

PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒
(Per 21 CFR 801.109)

OR

Over-the-Counter Use _____
(Optional Format 1-2-96)

Susan Walker
(Division Sign-Off)
Division of General, Restorative
and Neurological Devices

510(k) Number K 012806