

NOV 06 2001

K013316

October 2, 2001

Sickle-Chex  
510(k) Notification

**Names**

Trade Name	Sickle-Chex
Common Name:	Assayed Sickle Cell Control
Classification Name:	Sickle Cell Test, 21 CFR 864.7825
Classification Number:	Not Available

**Establishment Registration Number**  
1950302

**Classification**

FDA has assigned Sickle Cell Test to Class II, performance standards.

**Actions Taken to Comply with FDCA Section 154**

No government mandatory or industry voluntary performance standards exist for this type of device.

**Representative Labels and Labeling**

Proposed vial labels and other labeling attached:

- a. Package insert/ assay sheet
- b. Vial labels (Table I)
- c. Outer package labeling (Table I)

Predicate product package insert and assay sheet

**Physical description of Device**

Sickle-Chex consists of Control materials for verifying performance of sickle cell hemoglobin testing systems. The controls contain human red blood cells and preservative suspension media packaged in a polyethylene dropper bottle with dispensing tip. Product fill is 2.5 ml per vial.

**Intended Use**

Sickle-Chex consists of two controls. The Positive Control is a suspension of human red blood cells in a preservative media. This product is intended for use as a quality control material for sickling hemoglobin test kits. The Positive Control will produce a positive test for sickling hemoglobin using standard sickle cell hemoglobin solubility test reagents. In like manner, the Negative Control contains normal human red blood cells that will produce a negative test for sickling hemoglobin using standard sickle cell hemoglobin solubility test reagents.

Sickle-Chex may also be used as a control in hemoglobin electrophoresis methods.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Mr. Paul Kittelson  
Quality Assurance/Regulatory Affairs  
Streck Laboratories, Ins.<sup>®</sup>  
7002 South 109<sup>th</sup> Street  
LaVista, NE 68128

NOV 06 2001

Re: k013316  
Trade/Device Name: Sickie-Chex  
Regulation Number: 21 CFR 864.7825  
Regulation Name: Sickie cell test  
Regulatory Class: Class II  
Product Code: GHM  
Dated: October 3, 2001  
Received: October 4, 2001

Dear Mr. Kittelson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 -

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive, slightly slanted style.

Steven I. Gutman, M.D., M.B.A.  
Director  
Division of Clinical Laboratory-Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

NOV 06 2001

K013316

Page 1 of 1

510(k) Number:

K013316

Device Name:

Sickle-Chex

Indications For Use:

Sickle-Chex is intended to be used as a sickle cell control in testing for the presence of hemoglobin S in solubility tests and hemoglobin electrophoresis.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_

(Optional Format 1-2-96)

Date: \_\_\_\_\_

Joseph Bantua  
(Division Sign-Off)

Division of Clinical Laboratory Devices

510(k) Number K013316