

MAR 14 2002

EXHIBIT 2

SUMMARY OF SAFETY AND EFFECTIVENESS DATA

RESPONSE TO SMDA OF 1990

Kenneth J. Berk
80 Oakland Street
PO Box 780
Watertown, MA 02472 USA

Telephone: 617-926-6666
Fax: 617-926-6262

DEVICE

Trade Name: **PULPDENT Triple-Cure™ Reinforced Glass Ionomer Orthodontic Cement**
Classification Name: Adhesive, Bracket and Tooth Conditioner, Resin
FDA Product Code: DYH, 21 CFR Part 872.3750

PREDICATE DEVICES

GC Fuji Ortho LC (light-cure)
GC Fuji Ortho Self-cure
Pulpdent Band-Rite

DESCRIPTION AND INTENDED USE

PULPDENT Triple-Cure™ Reinforced Glass Ionomer Orthodontic Cement is a multi-purpose orthodontic cement that can be used as a bracket adhesive and as a band cement. It combines the advantages of glass ionomer cement with the added strength of dental resin technology. Triple-cure™ is a fluoride releasing, light-cure, self-cure, glass ionomer dental adhesive in powder and liquid form. It is used with an enamel etching gel to adhere orthodontic brackets to tooth surfaces.

COMPARISON WITH PREDICATE PRODUCTS:

PULPDENT Triple Cure Orthodontic Bracket Adhesive is substantially equivalent in composition and intended use to the predicate products listed above. Please see Exhibit 4 for the entire comparison.

SAFETY AND EFFECTIVENESS:

According to the NIH Technology Assessment Conference on *Effects and Side-Effects of Dental Restorative Materials*: "General usage of these materials over about 20 years indicates a high benefit-to-risk ratio...both composites and glass ionomers are relatively trouble-free. There is no evidence of short-term or long-term risk...There is no suspicion of any problems after virtually billions of procedures in the United States. In addition, the predicate products listed above have been given 510 (k) premarket approval as Class II Dental Devices under CFR 872.3750.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

MAR 14 2002

Mr. Kenneth J. Berk
Director
Puldent Corporation
80 Oakland Street
Watertown, Massachusetts 02472

Re: K014138

Trade/Device Name: Puldent Triple-Cure Reinforced Glass Ionomer Orthodontic
Regulation Number: 872.3750
Regulation Name: Adhesive, Bracket and Tooth Conditioner, Resin
Regulatory Class: II
Product Code: DYH
Dated: December 13, 2001
Received: December 17, 2001

Dear Mr. Berk:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.


You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4613. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Timothy A. Ulatowski
Director

Division of Dental, Infection Control
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

INDICATIONS FOR USE STATEMENT

510 (k) Number
(if known)

K014138

Device Name: **PULPDENT Triple-Cure™ Reinforced Glass Ionomer Orthodontic Cement**

Indications for Use:

PULPDENT Triple-Cure™ Reinforced Glass Ionomer Orthodontic Cement is a multi-purpose orthodontic cement that can be used as a bracket adhesive and as a band cement. It combines the advantages of glass ionomer cement with the added strength of dental resin technology. Triple-Cure™ is a fluoride releasing, light-cure, self-cure, glass ionomer dental adhesive in powder and liquid form. It is used with an enamel etching gel to adhere orthodontic brackets to tooth surfaces


Please do not write below this line. Continue on another page if needed.

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use _____
(Per 21 CFR 801.109)

or

Over-The-Counter Use



(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices
510(k) Number K014138