

JAN 15 2002

**510(K) SUMMARY**  
(as Required by 21 CFR § 807.92)

K014161

**A. Submitters Information**

**Submitter's Name:** St. Jude Medical, Inc  
Cardiac Surgery Division

**Address:** St. Jude Medical, Inc.  
One Lillehei Plaza  
St. Paul, MN 55117

**Contact Name** William McKelvey  
Regulatory Affairs Coordinator  
St. Jude Medical, Inc.  
Bus : (651) 415-7029  
Fax : (651) 766-3049  
Email: wmckelvey@sjm.com

**Submission Prepared:** December 17, 2001

**B. Device Information**

**Proprietary Name:** SJM® Tailor™ annuloplasty ring  
(Tailor ring)

**Common or Usual Name:** Flexible Annuloplasty Ring  
"C" ring

**Classification:** Pre-amendment Class II CFR § 870.3800  
Cardiovascular Prosthetic Devices,  
Annuloplasty Ring (revised April 10, 2001)

**Predicate Device:** St. Jude Medical considers The Tailor ring,  
model TARP to be substantially equivalent  
to the Tailor ring model TARN.

**Device Description** The Tailor ring is a flexible ring fabricated  
with a medical grade silicone rubber core,  
surrounded by polyester fabric. The full  
ring may be cut to form a partial or "C" ring.

**Intended Use:**

The Tailor ring is indicated for use in repair of diseased or damaged mitral or tricuspid heart valves that are determined by the physician to be repairable and do not require replacement.

**C. Comparison of Required Technological Characteristics**

SJM considers the Tailor ring, model TARP to be substantially equivalent in configuration, function and intended use to the Tailor ring, model TARN. The table below is a comparison of the equivalency characteristics between the two devices.

Characteristic	Equivalency
a. Product Labeling	Substantially Equivalent
b. Intended Use	Identical
c. Physical Characteristics	Different (Holder and Handle only)
d. Anatomical Sites	Identical
e. Target Population	Identical
f. Performance Testing	Substantially Equivalent
g. Safety Characteristics	Substantially Equivalent

**D. Summary of Non-Clinical Tests**

The testing for the Tailor ring model TARN (predicate) is included in the premarket notification (K000119). The following tests have been performed on the Tailor ring model TARP to insure substantial equivalence with the predicate.

**New Holder/Handle Configuration**

1. Physical Testing
  - Holder to Handle connection
  - Holder assembly
  - Ring assembly to holder
2. Microbiological Testing
  - Biocompatibility
  - Sterility Assurance
3. Manufacturing Process Validation



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JAN 15 2002

Mr. William McKelvey  
Regulatory Affairs Coordinator  
St. Jude Medical, Inc.  
One Lillehei Plaza  
St. Paul, MN 55117

Re: K014161  
Trade Name: SJM® Tailor™ Annuloplasty Ring, Model TARP-(size)  
Regulation Number: 21 CFR 870.3800  
Regulation Name: Annuloplasty Ring  
Regulatory Class: Class II (two)  
Product Code: KRH  
Dated: December 17, 2001  
Received: December 19, 2001

Dear Mr. McKelvey:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set


Page 2 - Mr. William McKelvey

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4586. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Bram D. Zuckerman, M.D.  
Acting Director  
Division of Cardiovascular  
and Respiratory Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number (if known): K014161


Device Name: SJM® Tailor™ annuloplasty ring

Indications for Use:

The SJM® Tailor™ annuloplasty ring is indicated for use in the repair of a mitral or tricuspid valve that is diseased or damaged due to acquired or congenital processes.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF ANOTHER PAGE IS NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
Division of Cardiovascular & Respiratory Devices  
510(k) Number: K014161

Prescription Use X

or

Over-The-Counter Use \_\_\_\_\_

Per 21 CFR 801.109)

Optional Format 1-2-96)