

K014207

Page 1 of 1

**510(k) SUMMARY OF SAFETY AND EFFECTIVENESS**  
in Accordance with SMDA of 1990**REUSABLE INSTRUMENT SYSTEM**

December 19, 2001

**COMPANY:** Aesculap®, Inc.  
3773 Corporate Parkway  
Center Valley, PA 18034

**CONTACT:** Lisa M. Millington, Regulatory Associate  
800-258-1946 (phone) x5072  
610-791-6882 (fax)  
[lisa.millington@aesculap.com](mailto:lisa.millington@aesculap.com) (email)

**TRADE NAME:** Reusable Instrument System

**COMMON NAME:** Disposable / Reusable Instrument System

**DEVICE CLASS:** Class II

**PRODUCT CODE:** GEI

**CLASSIFICATION:** 878.4800  
Device, Electrosurgical, Cutting & Coagulation & Accessories

**REVIEW PANEL:** General & Plastic Surgery

**INTENDED USE**

Aesculap's Reusable Instrument System is designed to cut, dissect, manipulate and/or cauterize various tissue during endoscopic/laparoscopic, general, vascular, gynecological, and thoracic surgical procedures.

**DEVICE DESCRIPTION**

Aesculap's Reusable Instrument System is comprised of non-sterile, reusable handles which are capable of monopolar electrocautery, via a underside connection terminal, reusable jaw inserts which come in different sizes and one single use endoscopic scissors. The jaw inserts include the blades, shaft and locking nut. The reusable scissors insert is supplied non-sterile, and the single use scissors insert is supplied sterile. The reusable and single use jaw inserts are intended to be use with either the non-ratcheted or ratcheted handle.

**PERFORMANCE DATA**

No performance standards have been promulgated under Section 514 of the Food, Drug and Cosmetic Act for these devices. The Reusable Instrument System, however, conform to the following electromedical standard: IEC/EN 60601-2-2 and IEC/EN 60601-2-18.

**SUBSTANTIAL EQUIVALENCE**

The Reusable Instrument System are substantially equivalent in their intended use, material composition, labeling, design and basic operating principles to the following predicate devices:

- Microline:
  - 3 MM Selec-Tip Laparoscopic Instrument System (K980758)
  - "Selecta-tip" Laparoscopic Instrument System (K970826)
- Weck's Endoscopic Instrument System (Disposable/Reusable) (K926203)



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

MAR 21 2002

Ms. Lisa M. Millington  
Regulatory Associate  
Aesculap, Inc.  
3773 Corporate Parkway  
Center Valley, Pennsylvania 18034

Re: K014207  
Trade/Device Name: Repposable Instrument System  
Regulation Number: 878.4400  
Regulation Name: Electrosurgical, Cutting and Coagulation Devices  
and Accessories  
Regulatory Class: II  
Product Code: GEI  
Dated: December 20, 2001  
Received: December 21, 2001

Dear Ms. Millington:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Lisa M. Millington

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

*Miriam C. Provost*  
for Celia M. Witten, Ph.D., M.D.  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## INDICATIONS FOR USE STATEMENT

510(k) Number (if known): K014207

Device Name: Reposable Instrument System

## Indication for Use:

Aesculap's Reposable Instrument System is designed to cut, dissect, manipulate and/or cauterize various tissue during endoscopic/laposcopic general, vascular, gynecological, and thoracic surgical procedures.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒ or Over-the-Counter Use ☐  
(per 21 CFR 801.109)

(Optional Format 3-10-98)

Miriam C. Provost  
(Division Sign-Off)  
Division of General, Restorative  
and Neurological Devices

510(k) Number K014207