

FEB 01 2002

510(k) Summary

as required by CFR section 807.92(c)

I. General Information

Date: January 2, 2002

Applicant: Progeny, Inc.
1807 Barclay Blvd.
Buffalo Grove, Ill.
60089

Contact Person: Alan Krema

Telephone: 847-850-3800 x785

Fax: 847-850-3800

II. Names

Device Name:

Trade Name: JB-70
Common Name: Intra Oral X-Ray System
Classification Name: 76 EHD – Unit, X-Ray, Extra oral with Timer

III. Predicate Devices

Gendex Gx-770
Gendex 765DC

IV. Product Description

The Progeny, Inc. JB-70 Intraoral Dental X-Ray System is an extraoral source of x-rays for intraoral radiographs in dental diagnostic radiography.

The Progeny, Inc. JB-70 Intraoral Dental X-Ray System consists of the following main components:

- X-ray tubehead
- Rotating yoke for tubehead mounting
- Articulation arm
- Horizontal extension arm
- Electronic control unit
- Wall mount
- 8 inch cone

Optional Components:

- 12 inch cone
- 8 ft. coil cord with exposure switch

The Power supply is regulated to provide a fixed 70 kVp at a fixed tube current of 8 mA. Predefined exposure times may be selected via the operator control panel. The range of exposure times is 0.05 to 1.65 seconds.

V. Indications for Use / Rationale for Substantial Equivalence

The JB-70 Intraoral Dental X-Ray System is to be used as an extraoral source of x-rays in Dental radiography.

The JB-70 shares the same indications for use, materials, design, operational and functional features and is therefore substantially equivalent to the predicate devices listed in section III of this summary.

There are several major independent manufacturers of Intra Oral diagnostic Radiographic Systems on the U.S. market. One is the Dentsply Gx770. The 510(k) number is K935046 . The classification of the Dentsply device is listed as product code 76EHD.

The other currently marketed device is the Gendex 765, manufactured by Dentsply, Inc. the 510(k) number is K992610. The classification of the Gendex 765 is listed as product code 76EHD .

Labeling for the currently marketed devices is included as Appendix B.

Comparison Table:

Characteristic	Gx 770	Gendex 765	Progeny JB-70
Fixed Kv	70 kVp	65 kVp	70 kVp
Fixed mA	7.0 mA	7.5 mA	8.0 mA
Operator Technique Selection	Exposure Duration Only	Exposure Duration Only	Exposure Duration Only
Suspension	Articulated Arm	Articulated Arm	Articulated Arm
X-Ray Field Size	6.5 cm diameter	6.0 cm diameter	6.5 cm diameter
Focal Distance	8 in. cone 12 in. cone	8 in. cone	8 in. cone 12 in. cone
Operator Exposure Control	Deadman Switch	Deadman Switch	Deadman Switch

VI. Safety and Effectiveness Information

Safety and effectiveness is demonstrated by:

- Performance testing and verification to meet product specifications.
- Software testing to validate software design and performance.
- Hazard analysis and risk level assessment.
- Same indications for use as predicate devices.

All of the above steps and evaluations combine to demonstrate that the JB-70 Intraoral Dental X-Ray System is safe and effective when the device is used as labelled.

VII. Conclusion

The Progeny, Inc. JB-70 Intraoral Dental X-Ray System is determined to be substantially equivalent to the predicate devices, the Gendex Gx-770, and the Gendex 765DC. The JB-70 shares the same indications for use, materials, design, operational and functional features to the currently marketed predicate devices listed in section III of this summary. The JB-70 Intraoral Dental X-Ray System is safe and effective when the device is used as labelled.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

FEB 01 2002

Mr. Alan Krema
Director of Product Development
Progeny, Inc.
1407 Barclay Blvd.
BUFFALO GROVE IL 60089

Re: K020070
Trade/Device Name: JB-70 Intraoral X-Ray System
Regulation Number: 21 CFR 872.1800
Regulation Name: Extraoral source x-ray system
Regulatory Class: II
Product Code: 90 EHD
Dated: January 2, 2002
Received: January 9, 2002

Dear Mr. Krema:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.


This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

8xx.1xxx	(301) 594-4591
876.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4616
884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx	(301) 594-4616
892.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4654
Other	(301) 594-4692

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Nancy C. Brogdon
Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

510(k) Number (if known): _____

Device Name: _____

Indications For Use:


Statement of Indications for Use:

Intended Use for the Progeny JB-70:

The intended use of the Progeny JB-70 Intra Oral X-Ray System is to act as a diagnostic source for radiographic dental imaging.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)
Division of Reproductive, Abdominal,
and Radiological Devices
510(k) Number 4020070

Prescription Use ☒
(Per 21 CFR 801.109)

OR

Over-The-Counter Use _____

(Optional Format 1-2-96)