

FEB 1 2 2002

K020295 1/2

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS***The NexGen LPS TM Monoblock Tibia***

Submitter Name: Implex Corp.

Submitter Address: 80 Commerce Drive
Allendale, New Jersey 07401-1600

Contact Person: Les Heimann

Phone Number: (201) 818-1800

Fax Number: (201) 995-9722

Date Prepared: January 24, 2001

Device Trade Name: The NexGen LPS Trabecular Metal Monoblock Tibia

Device Common Name: Tibial Components

Classification Number and Name: 21 CFR § 888.3560

**Substantial
Equivalence:**

The term "substantial equivalence" as used in this 510(k) notification is limited to the definition of substantial equivalence found in the Federal Food, Drug and Cosmetic Act, as amended and as applied under 21 CFR 807, Subpart E under which a device can be marketed without premarket approval or reclassification. A determination of substantial equivalency under this notification is not intended to have any bearing whatsoever on the resolution of patent infringement suits or any other patent matters. No statements related to, or in support of substantial equivalence herein shall be construed as an admission against interest under the US Patent Laws or their application by the courts.

Device Description:

The NexGen LPS TM Monoblock Tibia is manufactured from Trabecular Metal (Hedrocel Porous Tantalum) with direct compression molded ultra-high molecular weight polyethylene (UHMWPE).

These tibial components are intended for use with Zimmer NexGen LPS Femoral Components.

510(k) Summary (Continued)

Indications for Use: The NexGen LPS TM Monoblock Tibia is intended for use where severe degeneration, trauma, or other pathology of the knee joint indicates cemented total knee arthroplasty.

Conclusion: The NexGen LPS TM Monoblock Tibia is substantially equivalent to the identified predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

FEB 20 2002

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. Les Heimann
Director of Regulatory Affairs
Implex Corp.
80 Commerce Drive.
Allendale, New Jersey 07401

Re: K020295

Trade Name: NexGen LPS TM Monoblock Tibia

Regulation Number: 21 CFR 888.3560

Regulation Name: Knee joint patellofemorotibial polymer/metal/polymer semiconstrained
cemented prosthesis

Regulatory Class: II

Product Code: JWH

Dated: January 25, 2002

Received: January 28, 2002

Dear Mr. Heimann:

This letter corrects our substantially equivalent letter of February 12, 2002, regarding the trade name.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent [(for the indications for use stated in the enclosure)] to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,


for Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if
known):

K020295

Device Name:

The NexGen LPS TM Monoblock Tibia

Indications For Use:

The NexGen LPS Trabecular Metal Monoblock Tibia is intended for use where severe degeneration, trauma, or other pathology of the knee joint indicates cemented total knee arthroplasty.

Miriam C. Provost

(Division Sign-Off)

Division of General, Restorative
and Neurological Devices

510(k) Number K020295

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH; Office of Device Evaluation (ODE)

Prescription
Use

(Per 21 CFR 801.109)

OR...

Over-The-
Counter Use

(Optional Format 1-2-96)