

FEB 22 2002

K020436  
Page 1 of 2

EXHIBIT 2

510(k) Summary of Safety and Effectiveness

|   |  |
|---|--|
| SEDECAL SA<br>Pelaya 9- Poligono Industrial Rio De Janeiro<br>28110 -Algete<br>Madrid Spain<br>Tel (34) 91-628 0544/91-628 1592<br>Fax (34) 91-628 0574<br>(Foreign Manufacturer) | SEDECAL USA, Inc.<br>2910 N. Arlington Heights Rd.<br>Arlington Heights Illinois 60006<br>Tel 847-394-6960<br>Fax 847-394-6966<br>(Initial Importer)<br>Contact: Gary Fromberg |
|---|--|

February 5, 2002

510(k) Summary of Safety and Effectiveness

1. Identification of the Device:  
Proprietary-Trade Name: Models SP-HF-2.8 and SP-HF-4.0 Portable X-ray Units  
Classification Name: Mobile X-ray system, Product Code 90 IZL  
Common/Usual Name: Portable general purpose diagnostic X-ray Unit.
2. Equivalent legally marketed devices: This product is similar in function to the MinXray HF100H (a pre-amendments device)
3. Indications for Use (intended use) Models SP-HF-2.8 and SP-HF-4.0 are intended for use by a qualified/trained doctor or technician on both adult and pediatric subjects for taking diagnostic radiographic exposures of the skull, spinal column, chest, abdomen, extremities, and other body parts. Applications can be performed with the patient sitting, standing, or lying in the prone or supine position.
4. Description of the Device: Models SP-HF-2.8 and SP-HF-4.0 are a portable units which operate from 120 V 50-60~ AC. The unit utilizes a newly designed high frequency inverter and can be either mounted to a tripod or support arm or can be hand held. The usual safety precautions regarding the use of x-rays must be observed by the operator.
5. Safety and Effectiveness, comparison to predicate device. The results of bench and user testing indicates that the new device is as safe and effective as the predicate device.

**6. Substantial Equivalence Chart, Models SP-HF-2.8 and SP-HF-4.0**

| Characteristic                        | MinXray HF100H  | Models SP-HF-2.8   | Model SP-HF-4.0  |
|---------------------------------------|---|--|--|
| Intended Use:                         | Potable general purpose diagnostic X-ray unit   | SAME   | SAME   |
| Physical characteristics:             |   |  |  |
| Size                                  | 8.75"H x 9.5" W x 15.35" D.   | 8.7 H x 10.4" W x 16.5" D  | 8.7 H x 10.4" W x 16.5" D  |
| Weight                                | 40.9 lbs  | 33 lb  | 33 lb  |
| Energy Source:                        | 120 v 50-60~ AC   | 90 to 285 VAC (50-60 Hz)   | 90 to 285 VAC (50-60 Hz)   |
| User Interface                        | Up-Down pushbuttons for three kVp selections and exposure time selections with LED indicators | Up-Down pushbuttons for kVp and mAs. kVp adjustable in 1 kVp steps | Up-Down pushbuttons for kVp and mAs. kVp adjustable in 1 kVp steps |
| Exposure times                        | 0.08-2.00 Sec. In 192 steps   | 0.002-10 sec 38 steps  | 0.001-10 sec 41 steps  |
| Ma.                                   | 20 mA   | 5, 6.4, 8, 10, 12.5, 16, 20, 25, 32, 40, 50                        | 5, 6.4, 8, 10, 12.5, 16, 20, 25, 32, 40, 50, 64, 80, 100           |
| KvP                                   | 100 KvP   | 110 KvP  | 115 KvP  |
| Standards and Safety characteristics: |   |  |  |
| Performance Standard                  | 21 CFR 1020.30  | SAME   | SAME   |
| Electrical safety:                    | UL 2601, IEC 60601-1  | SAME   | SAME   |

**7. Conclusion**

After analyzing both bench and user testing data, it is the conclusion of Sedecal that the Models SP-HF-2.8 and SP-HF-4.0 are as safe and effective as the predicate device, have few technological differences, and has no new indications for use, thus rendering it substantially equivalent to the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

FEB 22 2002

SEDECAL USA, Inc.  
% Mr. Daniel Kamm, P.E.  
Kamm & Associates  
PO Box 7007  
DEERFIELD IL 60015

Re: K020436  
Trade/Device Name: Models SP-HF-2.8 and SP-HF-4.0  
Portable X-ray Unit  
Regulation Number: 21 CFR 892.1720  
Regulation Name: Mobil x-ray system  
Regulatory Class: II  
Product Code: 90 IZL  
Dated: February 6, 2002  
Received: February 8, 2002

Dear Mr. Kamm:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.


This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

|                                  |                |
|----------------------------------|----------------|
| 8xx.1xxx                         | (301) 594-4591 |
| 876.2xxx, 3xxx, 4xxx, 5xxx       | (301) 594-4616 |
| 884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx | (301) 594-4616 |
| 892.2xxx, 3xxx, 4xxx, 5xxx       | (301) 594-4654 |
| Other                            | (301) 594-4692 |

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

j) Indications for Use

510(k) Number K020436

Device Name: Models SP-HF-2.8 and SP-HF-4.0 Portable X-ray Units

Indications for Use: Models SP-HF-2.8 and SP-HF-4.0 Portable X-ray Units are intended for use by a qualified/trained physician or technician on both adult and pediatric subjects for taking diagnostic x-rays.

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓ OR Over the Counter Use \_\_\_\_\_  
(Per 21 CFR 801.109)

Nancy Chiodon  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number K020436