

MAR 26 2002



K020796

**GE Medical Systems**  
General Electric Company  
P O Box 414 Milwaukee, WI 53201

## **510(k) SUMMARY OF SAFETY AND EFFECTIVENESS**

This 510(k) summary of safety and effectiveness information is submitted in accordance with the requirements of 21 CFR Part 807.87(h).

### **Submitter**

Larry A. Kroger, Ph.D.  
Senior Regulatory Program Manager  
Telephone: (262) 544-3894, FAX: (262) 544-3863  
Date Prepared: February 8, 2002

### **PRODUCT IDENTIFICATION**

Name: CardIQ Analysis II

Classification Name: Accessory to Computed Tomography System

#### **Classification**

Panel: 892 - Radiology

#### **Classification**

Number: 892.1750

Manufacturer : General Electric Medical Systems  
283, rue de la Miniere  
78533 Buc Cedex, FRANCE

Distributor: General Electric Medical Systems, Milwaukee, WI

**Marketed Devices** CardIQ Analysis II is substantially equivalent to the device listed below:

Model: CardIQ Option 1.0  
Manufacturer: General Electric Medical Systems, Milwaukee, WI  
510(k) #: K003408

### **Device Description:**

CardIQ Analysis II is a post processing software option that can be used in the analysis of CT angiography images to display structures of the heart in a MIP, reformat or volume rendering view. The software has the ability to measure the diameter of the vessel, Hounsfield Units within coronary arteries, or plaque density within a vessel. Measurements, densities and images can all be printed to reports or saved to the AW workstation.

**Indications for Use:**

CardIQ Analysis II is a post processing software option for Advantage Windows Workstation (AW) Platform. This product can be used for the analysis of CT angiography images and for the assessment of plaque densities with the coronary arteries. It provides a number of display, measurement and batch filming features. The product can be used to aid trained physicians for visualizing and assessing cardiac anatomy, and coronary vessels.

**Comparison with Predicate:**

CardIQ Analysis II is a software option for CT Scanners. Features of this software package are substantially equivalent to the following device:

Device Name	FDA Clearance Number
CT CardIQ Option 1.0	K003408

**Adverse Effects on Health :**

The potential hazards are identified in a risk management summary (hazard analysis) and are controlled by:

- Software Development, Validation and Verification Process to ensure performance to specifications, Federal Regulations and user requirements.
- Adherence to industry and international standards.

**Conclusions:**

CardIQ Analysis II does not result in any new potential safety risks and performs as well as devices currently on the market. GE considers features of the CardIQ Analysis II to be equivalent to those of CT CardIQ Option 1.0 (K003408).



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

MAR 26 2002

General Electric Medical Systems  
% Mr. Wolfram Gmelin  
Technical Manager  
TÜV Rheinland of North America  
12 Commerce Road  
NEWTOWN CT 06470

Re: K020796  
Trade/Device Name: Card IQ Analysis II  
CT Scanner Post-Processing Software  
Regulation Number: 21 CFR 892.1750  
Regulation Name: Computed tomography x-ray system  
Regulatory Class: II  
Product Code: 90 JAK  
Dated: February 1, 2002  
Received: March 12, 2002

Dear Mr. Gmelin:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.


This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

8xx.1xxx	(301) 594-4591
876.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4616
884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx	(301) 594-4616
892.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4654
Other	(301) 594-4692

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## STATEMENT OF INTENDED USE

510(k) Number (if known): K020796

Device Name: CardIQ Analysis II

### Indications For Use:

CardIQ Analysis II is a post processing software option for Advantage Windows Workstation (AW) Platform. This product can be used for the analysis of CT angiography images and for the assessment of plaque densities with the coronary arteries. It provides a number of display, measurement and batch filming features. The product can be used to aid trained physicians for visualizing and assessing cardiac anatomy, and coronary vessels.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE  
IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒                       
(Per 21 CFR 801.109)

-OR-

Over-The-Counter Use                     

Nancy C. Brogdon  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number K020796