

APR - 4 2002



CORPORATE HEADQUARTERS

**SUMMARY OF SAFETY AND EFFECTIVENESS**

K020970  
page 1 of 2

**Sponsor:** Biomet, Inc.  
56 East Bell Drive  
P.O. Box 587  
Warsaw, IN 46581-0587

**Contact Person:** Tracy J. Bickel  
(574) 267-6639

**Proprietary Name:** Biomet® Oncology Salvage System

**Common Name:** OSS-K Diaphyseal Segments

**Classification Name:** Prosthesis, Hip, semi-constrained, metal/polymer,  
cemented (888.3350, 87JDI)  
Prosthesis, Knee, femorotibial, constrained, cemented,  
metal/polymer (888.3510, 87KRO)

**Substantially Equivalent Devices:** Oncology Salvage System (OSS)- K002757

**Device Description:** The Biomet® Oncology Salvage System offers a variety of component options for treatment of patients that require proximal femoral, distal femoral, total femur, or proximal tibial replacements, as well as, resurfacing components for the proximal tibial and distal femur.

There are twelve diaphyseal segment lengths available in 3cm, 4cm, and 5cm to 23 cm in 2cm increments. The diaphyseal segments are machined from wrought Ti-6Al-4V (ASTM F-1472). All the segments have a smooth satin finish.

**Intended Use:** The Biomet® Oncology Salvage System offers a variety of component options for treatment of patients that require proximal femoral, distal femoral, total femur, or proximal tibial replacements, as well as, resurfacing components for the proximal tibial and distal femur.

**Indications:**

- 1) Painful and disabled joint resulting from avascular necrosis, osteoarthritis, rheumatoid arthritis or traumatic arthritis
- 2) Correction of varus, valgus, post traumatic deformity
- 3) Correction of revision of unsuccessful osteotomy, arthrodesis, or previous joint replacement
- 4) Ligament deficiencies
- 5) Tumor resections
- 6) Treatment of non-unions, femoral neck, and trochanteric fracture of the proximal femur with head involvement, unmanageable using other techniques
- 7) Revision of previously failed total joint arthroplasty
- 8) Trauma

These devices are single use implants.

These devices are for cemented use only

MAILING ADDRESS  
P.O. Box 587  
Warsaw, IN 46581-0587

SHIPPING ADDRESS  
56 E. Bell Drive  
Warsaw, IN 46582

OFFICE  
219.267.6639

FAX  
219.267.8137

E-MAIL  
biomet@biomet.com

00064

Biomet, Inc.

Summary of Safety and Effectiveness – Page 2

*OSS-K Diaphyseal Segments*

K020970  
page 2 of 2

**Summary of Technologies:** The device's technological characteristics (materials, design, sizing, and indications) are similar to or identical to the predicate device.

**Non-Clinical Testing:** Surface finish analysis and mathematical calculations were performed to establish substantial equivalence.

**Clinical Testing:** None provided as a basis for substantial equivalence.

---

All trademarks are property of Biomet, Inc.

00065



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

APR - 4 2002

Ms. Tracy J. Bickel  
Regulatory Specialist  
Biomet Orthopedics, Inc.  
P.O. Box 587  
Warsaw, IN 46581-0587

Re: K020970

Trade/Device Name: Oncology Salvage System-K Diaphyseal Segments  
Regulation Number: 21 CFR 888.3350, 21 CFR 888.3510  
Regulation Name: Hip joint metal/polymer semi-constrained cemented prosthesis, Knee  
joint femorotibial metal/polymer constrained cemented prosthesis  
Regulatory Class: Class II  
Product Code: JDI, KRO  
Dated: March 25, 2002  
Received: March 26, 2002

Dear Ms. Bickel:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Tracy Bickel

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "for Mark A. Witten", is written over the typed name.

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative  
and Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

510(k) Number (if known): K020970

Device Name: ***OSS-k Diaphyseal Segments***

Indications for Use:

The Oncology Salvage System offers a variety of component options for treatment of patients that require proximal femoral, distal femoral, total femur, or proximal tibial replacements, as well as, resurfacing components for the proximal tibial and distal femur.

Indications:

- 1) Painful and disabled joint resulting from avascular necrosis, osteoarthritis, rheumatoid arthritis or traumatic arthritis
- 2) Correction of varus, valgus, post traumatic deformity
- 3) Correction of revision of unsuccessful osteotomy, arthrodesis, or previous joint replacement
- 4) Ligament deficiencies
- 5) Tumor resections
- 6) Treatment of non-unions, femoral neck, and trochanteric fracture of the proximal femur with head involvement, unmanageable using other techniques
- 7) Revision of previously failed total joint arthroplasty
- 8) Trauma

These devices are single use implants.

**These devices are for cemented use only.**

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF  
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use  
(Optional Format 1-2-96)

Mark J. Melanson  
(Division Sign-Off)  
Division of General, Restorative  
and Neurological Devices

00009

510(k) Number K020970