

MAY 15 2002

. 510(K) SUMMARY

Submitted by:

Organogenesis Inc.
150 Dan Road
Canton, Massachusetts 02021

Contact

Michael L. Sabolinski
Telephone: (781) 575-0775
Facsimile: (781) 575-0440

Date: April 25, 2002

Device:

Trade Name:	FortaPerm™
Common/Usual Name:	Surgical Mesh, Sling, Tissue Support Biomaterial
Classification Name:	Surgical Mesh (79FTM, 878.3300)
Regulatory Class:	Class II

Predicate Device:

The device is similar to predicate collagen-based surgical mesh and sling devices previously cleared for commercial distribution. The relevant predicate devices include GraftPatch (K970561) manufactured by Organogenesis Inc. and Surgisis Sling (K992159) manufactured by Cook Biotech, Incorporated.

Statement of Substantial Equivalence:

FortaPerm is substantially equivalent to the predicate devices, having similar intended use, technological characteristics, materials, physical construction and performance.

Intended Use:

FortaPerm is intended for implantation to reinforce and support soft tissue where weakness exists including, but not limited to the following procedures: pubourethral support, prolapse repair (urethral, vaginal, rectal and colon), reconstruction of the pelvic floor, bladder support, sacrocolposuspension, reconstructive procedures and tissue repair. By providing pubourethral support, FortaPerm may be used for the treatment of urinary incontinence resulting from urethral hypermobility or intrinsic sphincter deficiency. The device is intended for one-time use.

Device Description:

FortaPerm consists of a multi-laminate sheet predominantly of Type I porcine collagen. The device is supplied in sheet form in sizes ranging from 2 x 5 cm to 12 x 36 cm in sterile double laye packaging.

Performance Data:

FortaPerm was subjected to a panel of tests to assess biocompatibility, integrity, and performance. The device passed the requirements of all tests.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

MAY 15 2002

Mr. Michael L. Sabolinski, M.D.
Executive Vice President
Medical and Regulatory Affairs
Organogenesis, Inc.
150 Dan Road
Canton, MA 02021

Re: K021107
Trade/Device Name: FortaPerm™
Regulation Number: 878.3300
Regulation Name: Surgical Mesh
Regulatory Class: II
Product Code: FTM
Dated: April 4, 2002
Received: April 5, 2002

Dear Dr. Sabolinski:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Mr. Michael L. Sabolinski, M.D.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

Miriam C. Provost
for Celia Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

INDICATIONS FOR USE STATEMENT

Applicant: Organogenesis Inc.

510(k) Number (if known): K02 1107

Device Name: FortaPerm™

Indications For Use:

FortaPerm™ is intended for implantation to reinforce and support soft tissues where weakness exists, including but not limited to the following procedures: pubourethral support, prolapse repair (urethral, vaginal, rectal and colon), reconstruction of the pelvic floor, bladder support, sacrocolposuspension, reconstructive procedures and tissue repair. By providing pubourethral support, the FortaPerm™ may be used for the treatment of urinary incontinence resulting from urethral hypermobility or intrinsic sphincter deficiency.

The device is intended for one-time use.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE
IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Miriam C. Provost
(Division Sign-Off)
Division of General, Restorative
and Neurological Devices

Prescription Use X
(Per 21 CFR 801.109)

510(k) Number K02 1107 ~~OR~~ ~~Over-The-Counter Use~~ _____

(Optional Format 1-2-96)