



**U.S. FOOD & DRUG**  
ADMINISTRATION

November 18, 2019

Integra Lifesciences Corp.  
Diana Bordon  
Manager Of Special Projects  
311 Enterprise Dr.  
Plainsboro, New Jersey 08536

Re: K021792

Trade/Device Name: Bilayer Matrix Wound Dressing  
Regulatory Class: Unclassified  
Product Code: KGN  
Dated: May 30, 2002  
Received: May 31, 2002

Dear Diana Bordon:

This letter corrects our substantially equivalent letter of August 14, 2002.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

**Kimberly Ferlin -S**

Kimberly M. Ferlin, Ph.D.  
Assistant Director (acting)  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

**Bilayer Matrix Wound Dressing  
510(K) SUMMARY**

K021792

**Submitter's name and address:**

Integra LifeSciences Corporation  
311 Enterprise Drive  
Plainsboro, NJ 08536 USA

**AUG 14 2002**

**Contact person and telephone number:**

Diana M. Bordon  
Manager, Regulatory Affairs,  
(609) 275-0500

**Date:** May 23, 2002

**Name of the device:**

Proprietary Name: Bilayer Matrix Wound Dressing  
Common Name: Wound Dressing  
Classification Name: Dressing, Product Code 79FRO

**Substantial Equivalence:**

Bilayer Matrix Wound Dressing is substantially equivalent in function and intended use to the following products which have been cleared to market under Premarket Notifications 510(k): Oasis™ SIS Wound Dressing II (K993948), Fortaderm™ Wound Dressing (K014129), VitaChoice™ Wound Dressing (K896455) and Biobrane® II Temporary Wound Dressing (K896110).

**Intended Use:**

Bilayer Matrix Wound Dressing is indicated for the management of wounds including: partial and full thickness wounds, pressure ulcers, venous ulcers, diabetic ulcers, chronic vascular ulcers, surgical wounds (donor sites/grafts, post-Moh's surgery, post-laser surgery, podiatric, wound dehiscence), trauma wounds (abrasions, lacerations, second-degree burns, and skin tears) and draining wounds. The device is intended for one-time use.

**Device Description:**

Bilayer Matrix Wound Dressing is an advanced woundcare device comprised of a porous matrix of cross-linked bovine tendon collagen and glycosaminoglycan and a semi-permeable polysiloxane (silicone) layer. The semi-permeable silicone membrane controls water vapor loss, provides a flexible adherent covering for the wound surface and adds increased tear strength to the device. The collagen-glycosaminoglycan biodegradable matrix provides a scaffold for cellular invasion and capillary growth.

**Tests and Test Results**

Biocompatibility studies have demonstrated Bilayer Matrix Wound Dressing to be non-cytotoxic, non-pyrogenic, non-irritating, non-sensitizing, non-hemolytic and non-toxic.

**Conclusion**

Valid scientific evidence through biocompatibility and physical property testing provide reasonable assurance that Bilayer Matrix Wound Dressing is safe and effective under the proposed conditions of use, and is, with respect to intended use and technological characteristics, substantially equivalent to the predicate devices.

INDICATIONS FOR USE

Page 1 of 1

510(k) Number: K021792

Device Name: Bilayer Matrix Wound Dressing

Indications for Use:

Bilayer Matrix Wound Dressing is indicated for the management of wounds including: partial and full thickness wounds, pressure ulcers, venous ulcers, diabetic ulcers, chronic vascular ulcers, surgical wounds (donor sites/grafts, post-Moh's surgery, post-laser surgery, podiatric, wound dehiscence), trauma wounds (abrasions, lacerations, second-degree burns, and skin tears) and draining wounds. The device is intended for one-time use.

(PLEASE DO NOT WRITE BELOW THIS LINE—CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

[Signature]  
(Division Sign-Off)

Division of General, Restorative  
and Neurological Devices

Prescription Use X  
(Per 21 CFR 801.109)

Over-the-Counter Use \_\_\_\_\_

510(k) Number K021792 (Optional Format 1-2-96)