

K022037

	SPECIAL 510(k) SLEEP SYSTEM
SPECIAL 510(k) DEVICE MODIFICATION	JUNE 17, 2002 PAGE 49 of 53

Section F – 510(k) SUMMARY

SEP 10 2002

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of the SMDA 1990 and 21 CFR § 807.92.

Name: Cameron Mahon
Vice President, Customer Satisfaction

Address: **XLTEK**
2568 Bristol Circle
Oakville, Ontario
Canada, L6H 5S1

Telephone: (905) 829-5300

Fax: (905) 829-5304

E-mail: research@xltek.com

Common Names: Sleep Headbox

Classification Name: Electroencephalograph

Predicate Devices: XLTEK PSG-40 Polysomnography Headbox
K9991900

Masimo Radical Pulse Oximeter
K992340/K000126

XLTEK Ambulatory EEG
K982479

Description: The SLEEP Headbox is a digital polysomnograph headbox used in conjunction with XLTEK SLEEP Software to acquire and review sleep recordings (polysomnography).

Pg 1 of 2

 XLTEK ™ 	SPECIAL 510(k) SLEEP SYSTEM
SPECIAL 510(k) DEVICE MODIFICATION	JUNE 17, 2002 PAGE 50 of 53

Substantial Equivalence: The SLEEP Headbox is substantially equivalent to the XLTEK PSG-40 (K991900) incorporating Masimo pulse oximetry technology (K992340/K0000126). It comprises circuitry that is substantially equivalent to the XLTEK Ambulatory EEG (K982479). These modifications do not affect the predicate device's intended use, safety, or fundamental scientific technology.

Indications for Use:

The Sleep Headbox works in conjunction with Excel Tech Sleep software. This Sleep System is used to acquire and review sleep recordings (polysomnography) in research or clinical environments for:

- Digital recording of high-level output signals (such as EEG, respiratory and oximetry signals) from conventional polygraphic recorders, signal transducers or amplifiers.
- Selection of recorded signal sections for on-screen review, annotation and marking of sleep stages.
- Computer-assisted event marking and quantitative analysis of EEG, respiratory and oximetry signals.
- Computer-assisted reporting of simple measures obtained from the recorded signals (such as magnitude, time and frequency and simple statistical measures of marked events)

The SLEEP System is not intended to replace conventional devices or methods used for sleep monitoring in critical care or intraoperative settings.

The SLEEP System requires competent user input, and its output must be reviewed and interpreted by trained polysomnographers or trained medical professionals who will exercise professional judgment in using this information.

The SLEEP System does not make any judgment of normality or abnormality of the displayed signals or the result of an analysis. In no way are any of the functions represented as being in and of themselves diagnostic.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room - WO66-G609
Silver Spring, MD 20993-0002

Excel-Tech Ltd.
Sonja Markez
Regulatory Affairs
2568 Bristol Circle
Oakville, Ontario
Canada, L6H 5S1

APR - 9 2012

Re: K022037
Trade/Device Name: XLTEK Sleep System
Regulation Number: 21 CFR 882.1400
Regulation Name: Electroencephalograph
Regulatory Class: II
Product Code: OLV, GWQ
Dated (Date on orig SE ltr): August 9, 2002
Received (Date on orig SE ltr): August 12, 2002

Dear Ms. Markez:

This letter corrects our substantially equivalent letter of September 10, 2002.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Malvina B. Eydelman, M.D.
Director
Division of Ophthalmic, Neurological,
and Ear, Nose and Throat Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

	SPECIAL 510(k) SLEEP SYSTEM
SPECIAL 510(k) DEVICE MODIFICATION	JUNE 17, 2002 PAGE 51 of 53

Section G – INDICATIONS FOR USE

510(k) Number (if known): K022037

Device Name: SLEEP System

Indications for Use:

The Sleep Headbox works in conjunction with Excel Tech Sleep software. This Sleep System is used to acquire and review sleep recordings (polysomnography) in research or clinical environments for:

- Digital recording of high-level output signals (such as EEG, respiratory and oximetry signals) from conventional polygraphic recorders, signal transducers or amplifiers.
- Selection of recorded signal sections for on-screen review, annotation and marking of sleep stages.
- Computer-assisted event marking and quantitative analysis of EEG, respiratory and oximetry signals.
- Computer-assisted reporting of simple measures obtained from the recorded signals (such as magnitude, time and frequency and simple statistical measures of marked events)

The SLEEP System is not intended to replace conventional devices or methods used for sleep monitoring in critical care or intraoperative settings.

The SLEEP System requires competent user input, and its output must be reviewed and interpreted by trained polysomnographers or trained medical professionals who will exercise professional judgment in using this information.

The SLEEP System does not make any judgment of normality or abnormality of the displayed signals or the result of an analysis. In no way are any of the functions represented as being in and of themselves diagnostic.


(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21§ CFR 801.109)

OR Over-The Counter Use _____

(Optional Format 1-2-96)


 (Division Sign-Off)
 Division of General, Restorative
 and Neurological Devices

510(k) Number K022037