

OCT 04 2002

## 510(k) Summary

**Prepared:**

August 1, 2002

K023124

**Submitter:**

Company Name: Canon U.S.A., Inc. (U.S. designated agent for Canon Inc.)  
Company Address: One Canon Plaza  
Lake Success, NY 11042  
Contact Person: Sheila Driscoll, Senior Product Safety Engineer  
Phone Number: (516) 328-5602  
Fax number: (516) 328-5169

**Proposed Device:**

Reason For 510(k): New Model  
Manufacturer: Canon Inc.-  
Trade Name: Canon  
Model Name: CXDI-11 LANMIX MLT add. version  
Classification Name: 90MQB, Solid State X-ray Imager  
FDA 510(k)#: To be assigned

**Predicate Device:**

Manufacturer: Canon Inc.  
Trade Name: Canon  
Model Name: CXDI-11  
Classification Name: 90MQB, Solid State X-ray Imager  
FDA 510(k)#: K981556

**Description Of Device:**

The Canon X-ray digital camera model CXDI-11 LANMIX MLT add. version is used to directly capture and convert conventional projection X-ray images to digital images. A sub-sampled image can be displayed on a preview monitor for viewing. The diagnostic image can be transmitted through a DICOM compatible digital network for printing. The device provides digital image capture for conventional film/screen radiographic examinations.

The Canon X-ray Digital Camera CXDI-11 LANMIX MLT add. version is substantially equivalent to the Canon X-ray Digital Camera CXDI-11. It differs from the CXDI-11 in that the LAN MIX MLT software was added to make possible the Multi-objective Frequency Processing.

**Intended Use:**

Canon X-ray digital camera CXDI-11/ CXDI-11 LANMIX MLT add. version provide digital image capture for conventional film/screen radiographic examinations. the device is intended to replace radiographic film/screen systems in all general purpose diagnostic procedures.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room – WO66-G609  
Silver Spring, MD 20993-0002

Canon U.S.A., Inc.  
% Mr. Joseph Murnane  
Senior Staff Engineer  
Underwriters Laboratories, Inc.  
1285 Walt Whitman Road  
MELVILLE NY 11747-3081

AUG 23 2013

Re: K023124

Trade/Device Name: LANMIX MLT add. Version of Canon X-ray Digital Cameras  
Regulation Number: 21 CFR 892.1680  
Regulation Name: Stationary x-ray system  
Regulatory Class: II  
Product Code: MQB  
Dated: September 13, 2002  
Received: September 19, 2002

Dear Mr. Murnane:

This letter corrects our substantially equivalent letter of October 4, 2002.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into class II (Special Controls), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of

medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Parts 801 and 809), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely Yours,

A handwritten signature in black ink, appearing to read "Janine M. Morris", is written over the typed name and title.

Janine M. Morris  
Acting Director

Division of Radiological Devices  
Office of In Vitro Diagnostic Device  
Evaluation and Safety  
Center for Devices and Radiological Health

Enclosure

## Attachment 2

### Indications for Use Statement

510(k) Number  
(if known)

K023124

Device Name

LANMIX MLT add. version of Canon X-ray Digital Cameras

Indications for Use


LANMIX MLT add. version of Canon's X-ray Digital Cameras provides digital image capture for conventional film/screen radiographic examinations. The device is intended to replace radiographic film/screen systems in all general purpose diagnostic procedures.

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER  
PAGE IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒  
(Per 21 CFR 801.109)

OR Over-The-Counter Use ☐



(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K023124