

SEP 29 2003

510(k) Summary

K031225

Page 1/3

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR §807.92. Information supporting claims of substantial equivalence, as defined under the Federal Food, Drug and Cosmetic Act, respecting safety and effectiveness is outlined below.

Date Prepared:

April 12, 2003 (Revision 1 as of May 26, 2003)

Name of Submitter:

John E. Sawyer
Realistic Quality Solutions
C/O: GfE Medizintechnik GmbH
6254 South 1280 East
Salt Lake City, Utah 84121
(801) 685-7050

Realistic Quality Solutions has received authorization to submit this 510(k) on behalf of GfE Medizintechnik GmbH

Device Proprietary Name:

TIMESH also known as TIMESH-TC

Classification Name:

Surgical Mesh per 21 CFR 878.3300

Common/Usual Names:

Mesh, Metallic
Mesh, Non-Metallic
Polypropylene (Absorbable Mesh)
Polypropylene, (Non-Absorbable) Mesh

Predicate Device:

The subject device has the same indications for use as the predicates. The main differences are; changes in the construction of the mesh and the coating of the mesh.

K00112 2	Prolene Soft Mesh	Ethicon, Inc.
K90565 5	NonAbsorbable Polypropylene Surgical Mesh	United States Surgical Corp.
K91552 6	Coated SURGIPRO Polypropylene Surgical Mesh	United States Surgical Corp.

Device Description:

TiMESH, also known as TiMESH-TC, hereinafter known as TiMESH/TiMESH-TC" is a non-absorbable polypropylene mesh constructed of warp knitted fabric made from monofilament polypropylene fibers with all around proprietary coating. The TiMESH is constructed using a warp-knit (not woven, not knitted) process.

This construction along with the proprietary coating allows the TiMESH to have the strength, flexibility, as well as durability and surgical adaptability to allow for necessary tissue ingrowth.

The TiMESH with its current construction provides for elasticity in all directions. The construction of the mesh permits the use of any size mesh without it's unraveling. The wettability allows for the proper adaptation of the mesh to various stresses that can be encountered in the body. The product is sterilized by using ethylene oxide gas.

Intended Use:

The TiMESH/TiMESH-TC is specifically indicated for laparoscopic and open surgery for the repair of direct or indirect inguinal, femoral, umbilical, and incisional hernias; rectal, vaginal, and apical prolapses; and parietal reinforcement of tissues and abdominal wall repair. The TiMESH is specifically indicated for open surgery for the repair of direct or indirect inguinal and femoral hernias, parietal reinforcement of tissues, and abdominal wall repair. The TiMESH has the same indications as a combination of the predicate devices. TiMESH is a prescriptive device and should only be used by a licensed physician.

Indications Statement:

The TiMESH/TiMESH-TC The TiMESH, also known as TiMESH-TC, along with Ethicon Inc. Prolene Soft Mesh, U.S.S.C.'s NonAbsorbable Polypropylene Surgical Mesh and Coated SURGIPRO Polypropylene Surgical Mesh are intended to be used for the reinforcement of tissue during surgical repair. Thus, the TiMESH also known as TiMESH-TC mesh and all of the predicates have the same intended use. The TiMESH is specifically indicated for laparoscopic and open surgery for the repair of direct or indirect inguinal, femoral, umbilical, and incisional hernias; rectal, vaginal, and apical prolapses; and parietal reinforcement of tissues and abdominal wall repair.

510(k) Summary (Continued)**Indications Statement Cont'd:**

The TiMESH is specifically indicated for open surgery for the repair of direct or indirect inguinal and femoral hernias, pariental reinforcement of tissues, and abdominal wall repair. The TiMESH has the same indications as a combination of the predicate devices. TiMESH is a prescriptive device and should only be used by a licensed physician.

Technological Characteristics:

GfE Medizintechnik GmbH believes that the subject device is substantially equivalent to the predicate devices. The TiMESH/TiMESH-TC is constructed of the same materials as the predicate devices except for the coating. The Coating process allows the TiMESH to be more flexible than the predicates.

Performance Data:

The patient contact materials used in this device are common materials used in other medical devices as well as the predicate devices and have a long history of biocompatibility.

This device complies with the requirements of AAMI/ANSI/ISO 10993, Biological Evaluation of Medical Devices". In addition, bench testing was conducted in accordance with the FDA Guidance Document "Guidance for the Preparation of a Premarket Notification Application for a Surgical Mesh".

Conclusions:

Based on the information contained herein, we conclude that the changes are minor and that the subject device is substantially equivalent to the predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room - W066-G609
Silver Spring, MD 20993-0002

GfE Medizintechnik GmbH
% Mr. John Sawyer
President
Realistic Quality Solutions
6254 South 1280 East
Salt Lake City, Utah 84121

JUN 20 2012

Re: K031225

Trade/Device Name: TiMESH also known as TiMESH-TC
Regulation Number: 21 CFR 878.3300
Regulation Name: Surgical mesh
Regulatory Class: II
Product Code: FTL, OXJ
Dated: August 8, 2003
Received: August 12, 2003

Dear Mr. Sawyer:

This letter corrects our substantially equivalent letter of September 29, 2003

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set


Page 2 - Mr. John Sawyer

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,


Mark N. Melkerson
Director
Division of Surgical, Orthopedic
and Restorative Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k): K031225


Device Name: TiMESH

Indications for Use: The TiMESH is intended for the reinforcement of tissue during surgical repair. Specifically, the TiMESH is indicated for laparoscopic and open surgery for the repair of direct inguinal hernias, indirect inguinal hernias, femoral hernias, umbilical hernias, incisional hernias, parietal reinforcement of tissues and abdominal wall repair.

Prescription Use X AND/OR Over the Counter Use
(Per 21 CFR 801 Subpart D) (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)
Division of Surgical, Orthopedic,
and Restorative Devices

510(k) Number K031225
